Scrubber (Peeling ultrasónico)

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F335A

Scrubber (Ultrasonic peeling)

Treatment designed to deeply cleanse the skin, eliminating dead cells and horny layer; by means of a spatula that generates ultrasonic waves, and produces vibration that in contact with the tissue helps us to drag all the residues, as well as favouring the stimulation of the skin, excellent results in the elimination of blackheads, blackheads, helps to eliminate comedones, dirt, residues of cosmetics or make-up, grease and acne.

Scrubber F335A (Ultrasonic Peeling) has two ultrasonic peeling handpieces, a large spatula for the face and a smaller one to work around the eyes and small areas where we can not reach with the large one, and a bracelet for driving, to be placed by the customer.

Place the bracelet on the client's wrist, adjusting it and applying tonic or gel for the best conduction of the current, in this way you can start working. The spatula should be inclined at a 45° angle, if the spatula loses contact with the skin it will automatically stop working.

Choose the type of wave in STATE

Continuous: this is a wave that is not interrupted and is used for large surfaces, as in this way the wave emission is continuously emitted, which generates greater intensity of treatment, used for thicker skins, with a greater amount of impurities.

Pulsed: used in small areas, with less thickness and little accumulation of dirt as the wave emission is smaller, therefore, the effect will be less powerful and consequently the work will be less. Also used for areas where the bone tissue is more prominent.

Combined: it is an intermediate point between the two wave emissions.

Choose the treatment time, regulate the intensity in energy, always from less to more to the client's tolerance. The treatment time is approximately 10-15 minutes.

How to use the Scrubber (Ultrasonic Peeling)

- Connect the equipment correctly.
- Switch it on, and then the machine will be in 'stand by' mode, ready for use.
- Select the spatula you are going to work with, the type of wave, the desired "power" and "time" level. Please note (the standard default settings are "continuous wave", "maximum power output" and "15 minute timer"), however, you can change any of these settings during the course of the treatment phase.
- Press the "START / PAUSE" button, hold the handpiece between thumb and forefinger respectively. The circuit will close once the wristband is on the client and the spatula is in contact with the skin. The display will also show that the handpiece is active. Hold the spatula at a 45° downward angle and move it only in an upward direction to cleanse the skin, some tonic or gel should be used to work with the ultrasonic peel.
- Once you have finished, remove the handpiece from the skin and the wristband.
- You can also use the back of the spatula to rehydrate or soothe the skin by applying the product for penetration with ultrasound.

F335A

CAUTION

- Do not use this type of treatment on clients with the following health problems: heart disease, infectious disease, skin infection, skin allergies, clients with pacemakers or metal implants in their body. Note: All jewellery must be removed before performing this treatment.
- Do not allow the probe to be "on" for long periods of time without being used, as "overheating" may occur and damage the probe blade. 3. Do not rub the probe blade or use any abrasive device (sandpaper, metal wool, etc.) as this will damage the blade and possibly cause abrasions to the client's skin. Do not use the probe on any surface other than dermal tissue.
- Never allow the skin surface to dry out during treatment, add more product as needed.
- Take care not to damage the probe cable during use, this can be achieved by not stretching the cable too much, or ensuring that the cable is not kinked.
- During use, ensure that there is no excess product in the machine or main tower unit to avoid accidental damage. Never use corrosive substances with the ultrasonic peeling probe.
- If the unit does not operate as expected, or if you experience unusual operating circumstances, contact your sales agent immediately for further advice.
- Never attempt to open the machine or its accessories for inspection, as they do not contain "user serviceable parts". Inspection or repair by unauthorised persons is prohibited and will invalidate the manufacturer's warranty.
- If the power cord of the machine is damaged, it must be replaced with a "genuine" part from the manufacturer. Never operate machines with a damaged power cord. Contact your sales agent for more information.

F335A

Attention:

- The system operates with a positive ground circuit, for safety reasons. If the handpiece cannot detect the customer's skin (find a return circuit to the machine), the system will stop output, which will be indicated by a "flashing" spatula symbol on the LCD display. This may be caused by a bad connection, or due to a lack of conductive product on the customer's skin. Check that there are no loose connections, that sufficient product has been applied and restart (turn off and on) the machine. Once a good circuit has been re-established, the machine will restart and continue until the selected time has elapsed.
- If the above checks do not rectify the problem when the machine restarts, and the spatula touches the skin, there may be possible damage to the handpiece blade, or damage to the ultrasonic generators or their filaments. The ultrasonic generators and connecting filaments are sensitive parts of the equipment and should be treated as such. Damage may occur if the handpiece is dropped or hit with sufficient force, and is not "user serviceable" and, if you suspect damage, you should contact your sales agent immediately and seek further advice. Use only "genuine" replacement parts from your dealer. Failure to do so may damage the equipment or cause injury to the user. It will also invalidate the manufacturer's warranty.
- During ultrasonic peeling, the use of "essential oils" or "oil-based" products is prohibited. The skin must be kept moist at all times to avoid overheating of the skin surface. The best way is to keep the handpiece always moving while in contact with the skin. Never use high viscosity skin care products (especially masks) with any ultrasonic peel, as some penetration of the product will occur.
- Before carrying out an ultrasonic peeling treatment, it is advisable to "preset" the power (energy) level of the probe before use. This is so that the level is not too high and causes discomfort to the client, and to prevent the product from spreading to untreated areas.

Contraindications

- Pregnancy and breastfeeding.
- Very sensitive skin.
- Rosacea.
- Seborrhoeic dermatitis.
- Recent lesions, skin with wounds, infections, psoriasis, dermatitis or unhealthy skin.
- Shaved or depilated skin.
- Skin that has recently undergone deep chemical peeling, burns or surgery.
- Alterations in skin sensitivity, unable to differentiate between cold and heat.
- Poor healing.
- Previous history of keloids.
- Areas with tattoos on the face or permanent make-up.
- Skin conditions.
- Infectious processes or fungus.
- Severe acne.
- Caution when working on tissues with injected fillers, and in the case of clients who have recently had hyaluronic acid, Botox or any other type of product injected, you should wait a reasonable amount of time before having treatment.

Yubitzay Avila
Customer Training
T. +34 93 824 9020
y.avila@weelko.com
S. Yubitzay Avila

Weelko Barcelona S. L.
Pol. Ind. La Valldan
Camí de Garreta, 92
08600 Berga
Barcelona – Spain
T. +34 93 824 9020
info@weelko.com
www.weelko.com